## These are some helpful hints for safe use and disposal of this container.

## DO:

- 1. Do only use this container to dispose of your needles and lancets.
- 2. Do make sure your container is assembled correctly. Press firmly all around the lid (all corners if it is square) until you don't hear any more clicks. Test the lid to see if you can take it off.
- 3. Do keep the container away from children and pets. Children find the colour attractive and may be injured by the contents.
- 4. Do use the temporary closure on the lid to help prevent spilling when you are not using the bin.
- 5. Do make sure you have another container available before your current one is full.
- 6. Do when the container is full, close the lid until you hear it click. This means that it is locked.
- 7. Do fill in your name and the date when you lock the bin.
- 8. Do when the container is full and locked; take it to your nearest Hillingdon Community Health clinic or health centre.

## DON'T:

- 1. Don't put your needles and lancets in any other container.
- 2. Don't put any other waste in this container; it is only designed to take needles and lancets.
- 3. Don't put the container into your normal rubbish.
- 4. Don't leave the container anywhere that it could be easily knocked over.
- 5. Don't fill the container past the fill line. If you overfill the container, needles may begin to stick through the plastic wall and cause an injury to someone.